

MICROPREFS® MICROPLASTICS REFERENCE MATERIALS

ENVIRONMENTAL

Setting the Standard in Microplastics Testing with Chiron MicroPRefs®

Microplastics are plastic particles smaller than 5 millimeters that have become widespread environmental pollutants. These non-biodegradable particles accumulate in ecosystems and can enter the food and beverage supply, making them harmful to both the environment and human health.

To support accurate analysis, our scientists developed MicroPRefs® microplastics reference materials. Designed to detect the six most abundant plastics, they are available in two formats: a novel, easy-to-use soda tablet, and neat powder.

Chiron MicroPRefs® Features:

- Available polymers: PE, PP, PS, PET, PVC, and PC
- Particle size: 50-300 µm with irregular shapes that reflect real-world environmental samples
- Soda tablets: Ideal for particle counting and polymer identification by Optical Microscope, µ-FTIR, µ-Raman, LDIR
- Neat (powder) form: Optimized for polymer identification and purity analysis by Pyr-GC/MS or TED-GC/MS



gozeptometrix.com/4loz0a5

+1.800.274.5487

ANALYTICAL.CS@ZEPTOMETRIX.COM

ZeptoMetrix®
an antylia scientific company

Tablet

Part No.	Chiron Part No.	Description	No. of Particles	Particle Size	Unit of Measure	Units/ Pack	Price
24121-65	15459.X-TAB	MicroPRefs® Blank tablet			EA	1	EUR 117 USD 136 GBP 104
24121-65X3	15459.X-TABx3	MicroPRefs® Blank tablets (3/PK)			PK	3	EUR 294 USD 340 GBP 260
24121-68	15246.X-50/300-TAB	MicroPRefs® Polyethylene tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-68X3	15246.X-50/300-TABx3	MicroPRefs® Polyethylene tablets (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-71	15247.X-50/300-TAB	MicroPRefs® Polyethylene terephthalate tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-71X3	15247.X-50/300-TABx3	MicroPRefs® Polyethylene terephthalate tablets(3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-74	15248.X-50/300-TAB	MicroPRefs® Polystyrene tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-74X3	15248.X-50/300-TABx3	MicroPRefs® Polystyrene tablets (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-77	15249.X-50/300-TAB	MicroPRefs® Polypropylene tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-77X3	15249.X-50/300-TABx3	MicroPRefs® Polypropylene tablets (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-80	15250.X-50/300-TAB	MicroPRefs® Polyvinylchloride tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-80X3	15250.X-50/300-TABx3	MicroPRefs® Polyvinylchloride tablets (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-83	15251.X-50/300-TAB	MicroPRefs® Polycarbonate tablet	50-70 particles	50-300 µm	EA	1	EUR 214 USD 248 GBP 189
24121-83X3	15251.X-50/300-TABx3	MicroPRefs® Polycarbonate tablets (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 536 USD 620 GBP 473
24121-86	S-5401-50/300-TAB	MicroPRefs® Tablet: Mix of PE, PET, PS	70-100 particles	50-300 µm	EA	1	EUR 284 USD 329 GBP 251
24121-86X3	S-5401-50/300-TABx3	MicroPRefs® Tablets: Mix of PE, PET, PS (3/PK)	50-70 particles	50-300 µm	PK	3	EUR 712 USD 823 GBP 629

Tablet

Part No.	Chiron Part No.	Description	No. of Particles	Particle Size	Unit of Measure	Units/ Pack	Price
24121-89	S-5402-50/300-TAB	MicroPRefs® Tablet: Mix of PP, PVC, PC	70-100 particles	50-300 µm	EA	1	EUR 284 USD 329 GBP 251
24121-89X3	S-5402-50/300-TABx3	MicroPRefs® Tablets: Mix of PP, PVC, PC (3/ PK)	50-70 particles	50-300 µm	PK	3	EUR 712 USD 823 GBP 629

Neat

Part No.	Chiron Part No.	Description	Particle Size	Units/ Pack	Volume	Price
24121-92	15246.X-1MG	MicroPRefs® Polyethylene	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-92X10	15246.X-10MG	MicroPRefs® Polyethylene	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73
24121-93	15247.X-1MG	MicroPRefs® Polyethylene terephthalate	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-93X10	15247.X-10MG	MicroPRefs® Polyethylene terephthalate	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73
24121-94	15248.X-1MG	MicroPRefs® Polystyrene	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-94X10	15248.X-10MG	MicroPRefs® Polystyrene	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73
24121-95	15249.X-1MG	MicroPRefs® Polypropylene	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-95X10	15249.X-10MG	MicroPRefs® Polypropylene	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73
24121-96	15250.X-1MG	MicroPRefs® Polyvinylchloride	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-96X10	15250.X-10MG	MicroPRefs® Polyvinylchloride	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73
24121-97	15251.X-1MG	MicroPRefs® Polycarbonate	50-300 µm	1	1 mg	EUR 21 USD 25 GBP 19
24121-97X10	15251.X-10MG	MicroPRefs® Polycarbonate	50-300 µm	1	10 mg	EUR 82 USD 95 GBP 73

Neat

Part No.	Chiron Part No.	Description	Particle Size	Units/ Pack	Volume	Price
24121-98	15558.6-KIT-1MG	MicroPRefs® PE, PET, PS, PP, PVC, & PC Kit	50-300 µm	6	6 mg	EUR 135 USD 156 GBP 120
24121-98X10	15558.6-KIT-10MG	MicroPRefs® PE, PET, PS, PP, PVC, & PC Kit	50-300 µm	6	60 mg	EUR 472 USD 546 GBP 417

Prices are valid from September 1 to December 31, 2025



gozeptometrix.com/4loz0a5

+1.800.274.5487
ANALYTICAL.CS@ZEPTOMETRIX.COM

SPEX® | NSI | HPS | CHIRON

Analytical Reference Materials & PT by ZeptoMetrix®