



**NEW PRODUCT****7H-PFHpA**

In 2025, the European Committee for Standardization (CEN) approved EN 17681-1:2025 for the analysis of PFAS in textiles and textile products. Included in the updated list of target compounds is 7H-perfluoroheptanoic acid (**7H-PFHpA** or HPFHpA). Additionally, **7H-PFHpA** is listed as an analyte in the ISO/DIS 25652:2025 method for the analysis of PFAS in sediment, soil, sludge and waste. In response to these new methods, **Wellington** is pleased to announce the release of a certified reference standard for **7H-PFHpA**.

**7H-PFHpA**

Product Code	Description (methanol solution)	Units	Conc
NEW 7H-PFHpA	7H-Perfluoroheptanoic acid	1.2 mL	50.0 µg/mL

Full scan chromatogram, extracted mass spectrum, and selected MRM transitions of **7H-PFHpA**.**DISTRIBUTED BY:****BCP INSTRUMENTS**

12 avenue des Saules

69600 OULLINS, France

tel +33(0)4 72 49 72 65

www.bcp-instruments.com

contact@bcp-instruments.com

